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## Copyright

## Translations

Please read this manual carefully before using the Product, so as to protect **"the Technical Service Personnel"** and **"the Operator"** from any injury.

This appliance is a Class I medical device pursuant to European Directive on medical devices (MDD) 93/42/EEC (Annex IX) as amended and integrated.

The manufacturer declares that this Product complies with Annex I (Essential requirements) of Directive 93/42/EEC as amended and integrated and certifies such conformity by affixing the CE marking.

This installation manual is valid for the following models:

- |            |            |
|------------|------------|
| • LC101LRA | • LC104LRD |
| • LC102LRA | • LC106LRD |
| • LC103LRA | • LC107LRD |
| • LC104LRA | • LC109LRD |
| • LC105LRA | • LC110LRD |
| • LC106LRA | • LC111LRD |
| • LC107LRA | • LC112LRD |
| • LC109LRA | • LC113LRD |
| • LC101LRD | • LC114LRD |
| • LC102LRD | • LC115LRD |
| • LC103LRD | • LC116LRD |

The customer service is at your disposal in case of Product details, information concerning its use, identification of spare parts being required and for any other queries you might have concerning the appliance, for ordering spares and for matters relating to assistance and warranty.

- TECNO-GAZ S.p.A.
- Strada Cavalli, 4
- I-43038 Sala Baganza PR
- Tel.: ++39 0521 8380
- Fax: ++39 0521 833391
- E-mail: [info@tecnogaz.com](mailto:info@tecnogaz.com)

The reproduction and translation, including partial, of any part of this manual is forbidden without the written permission of TECNO-GAZ.

The original language of this manual is ITALIAN. For all translations, reference must be made to the original manual language.

## Index of contents

|                                                                                          |           |
|------------------------------------------------------------------------------------------|-----------|
| <b>KEY.....</b>                                                                          | <b>3</b>  |
| <b>1 GENERAL SAFETY INFORMATION.....</b>                                                 | <b>4</b>  |
| <b>2 General information .....</b>                                                       | <b>5</b>  |
| 2.1 Operator qualifications.....                                                         | 5         |
| 2.2 Packaging, transport, storage and characteristics of installation premises .....     | 5         |
| 2.3 Graphic signs and symbols used in the installation manual.....                       | 6         |
| 2.4 Graphic symbols used on packaging .....                                              | 6         |
| 2.5 Graphic symbols used on the Product.....                                             | 7         |
| 2.6 Warranty and liabilities.....                                                        | 8         |
| 2.7 Structural changes or variations.....                                                | 8         |
| <b>3 Instructions on how to prepare the premises mechanically and electrically .....</b> | <b>8</b>  |
| 3.1 Preparing the premises mechanically (ceiling and wall version).....                  | 8         |
| 3.2 Correctly wiring up the premises.....                                                | 9         |
| <b>4 Product installation .....</b>                                                      | <b>9</b>  |
| 4.1 Parts included in the package.....                                                   | 10        |
| 4.2 Ceiling and wall drilling instructions .....                                         | 10        |
| 4.3 Instructions for ceiling version of Product .....                                    | 12        |
| 4.3.1 Installation of the ceiling plate, bar, power supply .....                         | 12        |
| 4.3.2 Installation of structure to bar .....                                             | 14        |
| 4.3.3 Installation of swinging arm .....                                                 | 14        |
| 4.3.4 Installation of cupola.....                                                        | 16        |
| 4.3.5 Electrical connection .....                                                        | 18        |
| 4.3.6 Installation of ceiling cover .....                                                | 19        |
| 4.4 Installation of Product in mobile version.....                                       | 19        |
| 4.4.1 Installation of lamp stem.....                                                     | 19        |
| 4.4.2 Installation of swinging arm .....                                                 | 20        |
| 4.4.3 Installation of cupola.....                                                        | 21        |
| 4.4.4 Electrical connection .....                                                        | 21        |
| 4.5 Installation of Product in wall version.....                                         | 21        |
| 4.5.1 Installation of plate with wall switchboard .....                                  | 21        |
| 4.5.2 Installation of structure to plate.....                                            | 22        |
| 4.5.3 Installation of swinging arm .....                                                 | 24        |
| 4.5.4 Installation of cupola.....                                                        | 25        |
| 4.5.5 Electrical connection .....                                                        | 25        |
| 4.6 Protection fuses.....                                                                | 26        |
| 4.7 Handpiece fitting.....                                                               | 26        |
| 4.8 Mechanical adjustments .....                                                         | 26        |
| 4.9 First switch-on.....                                                                 | 26        |
| 4.10 Check the result of Product installation and testing before use .....               | 27        |
| <b>5 Troubleshooting .....</b>                                                           | <b>28</b> |

**PRODUCT****KEY**

The EM (Electro-Medical) EQUIPMENT to which this manual refers is a **SECONDARY SURGICAL LAMP (TREATMENT LAMP)**. For ease of description, in this manual this EM EQUIPMENT will be called "**Product**".

**OPERATOR**

Professional medical personnel (e.g., professional health personnel, expert person assisting the patient).

**RESPONSIBLE  
ORGANIZATION**

Entity accountable for the use and maintenance of an EM equipment or EM system (e.g., a hospital, an individual doctor or a non-expert person). Preparation and awareness are included in use.

**TECHNICAL  
SERVICE  
PERSONNEL**

The personnel (individuals or entity accountable to the responsible organization) that installs, assembles, maintains or repairs the equipment. Under certain circumstances, the safety of such persons depends on their knowledge and awareness and ability to take appropriate precautions when gaining access to hazardous parts partially. By way of example only, the following professional figures are deemed as SERVICE PERSONNEL:

- ⇒ Construction Engineer, Draughtsman, Building firm duly registered in the professional Register (for the masonry works)
- ⇒ Electrical Engineer Electro-technical expert qualified to work as an electrician (for the electrical works)

For the installation phase, as regards assembly operations only, a qualified person is deemed whosoever has attended a course organized by TECNO-GAZ or, alternatively, whosoever has carefully read the manual.

## **1 GENERAL SAFETY INFORMATION**

This manual is an integral part of the Product as indicated by European Directive 93/42/EEC and subsequent amendments and supplements. Read and keep this manual close to the Product.

- The Product is not suitable for use in explosion-risk areas.
- The Product is not suitable for use wherever there are inflammable mixes of anaesthetics with air, oxygen or N<sub>2</sub>O (laughing gas).
- The Product is not suitable for use in environments rich in oxygen and use is not intended in the presence of inflammable agents.

TECNO-GAZ disclaims all liability for any injury to persons or damage to things caused by the Product having been installed by persons who are not **"TECHNICAL SERVICE PERSONNEL"**.

The RESPONSIBLE ORGANIZATION is entirely responsible for Product installation activities; no costs or responsibilities relating to the installation and/or commissioning of the Product may therefore be traced back and/or in any case attributed to TECNO-GAZ.

The ceiling masonry works for Products to be installed on ceilings, and the electrical works for supplying power to the Product shall be carried out in a workmanlike manner by TECHNICAL SERVICE PERSONNEL to ensure these are sturdy and safe.

The electrical system in the premises must conform to IEC:60364-7-710 standard and any national regulations. A master switch must be installed with fuse or thermal magnetic circuit breaker to be able to interrupt power to the Product.

Installation  
Use  
Cleaning  
Routine maintenance  
Special maintenance

Assistance  
Demolition

## Packing

## Transport

## Storage

## Installation premises

## 2 General information

### 2.1 Operator qualifications

Qualification of personnel in charge of operating on the Product:  
Installer and/or qualified technician.  
Professional medical personnel.  
Properly trained medical and paramedical personnel.  
Qualified technician with required technical-professional skills.  
TECNO-GAZ or technical service personnel, the latter only for the fuse change.  
TECNO-GAZ or authorized Dealer.  
Comply with applicable laws on waste disposal. This product must not be disposed of in standard waste disposal bins. To avoid risks for the environment and health deriving from the dispersion of polluting substances in the environment, separate the various internal component parts such as iron, aluminium, plastic and electrical material, and dispose of these through authorized channels so as to ensure correct recycling.

### 2.2 Packaging, transport, storage and characteristics of installation premises

Cardboard boxes containing Product. Dispose of these in compliance with national directives applicable for waste disposal.

Product transport is done by land, sea or air according to the following characteristics:  
Temperature (°C): -15 / +60  
Humidity: 10 / 75 %  
Atmospheric pressure (h/Pa): 500 / 1060

The packaged Product must be stored (warehoused) in dry premises having the following characteristics:  
Temperature (°C): -15 / +60  
Humidity: 10 / 75 %  
Atmospheric pressure (h/Pa): 500 / 1060

The premises where the Product is started up must have the following characteristics:  
Temperature (°C): +10 / +40  
Humidity: 30 / 75 %  
Atmospheric pressure (h/Pa): 700 / 1060

## 2.3 Graphic signs and symbols used in the installation manual

The following safety measures must be put in place during Product installation, use and servicing.

To emphasize their importance, a number of safety precautions are repeated throughout the manual.

Follow the safety precautions before using or repairing the Product. Carefully abiding by the safety precautions improves the ability to use the Product safely and correctly and helps prevent incorrect maintenance which could be hazardous and cause damage. The safety measures are approximate and not exhaustive; the Operator, the Responsible Organization and the Technical Service Personnel must develop their capacities to upgrade and integrate them.

General warning signal

General mandatory code of conduct signal

General prohibition signal

## 2.4 Graphic symbols used on packaging

List of symbols on packaging boxes:

This side upwards

Fragile

Protect from rain

Do not stack packaging

Weight of packaging

Humidity to be complied with (max limit at top right and min limit at bottom left)

Pressure to be complied with (max limit at top right and min limit at bottom left)

Limit temperature (max limit at top right and min limit at bottom left)

## 2.5 Graphic symbols used on the Product

Below are the symbols to be found on the Product:

CE marking indicating the Product conforms to directive 93/42/EEC and subsequent amendments and supplements

Date of manufacture (month and year)

Manufacturer's address

Fuses used in the device

Comply with the instructions for use

Model

Serial number

Disposal

Protection earth

Neutral lead connection point

Line lead connection point

ON


OFF

Standby and switch-on

Pushing, resting on or lying on the product is forbidden

No stepping on surface

Only move the product after lowering the arm

CAUTION

## 2.6 Warranty and liabilities

TECNO-GAZ disclaims all liability as regards unreliable Product operation in the following cases:

- Installation, authorized modifications and repairs have not been performed by TECHNICAL SERVICE PERSONNEL.
- The Product has not been used for its intended purpose and in conformity with the operating instructions (see operation manual).
- The premises have not been approved for healthcare activities.
- The premises are not built in conformity with the law and applicable regulations.
- The electrical system in the premises is not in compliance with appropriate requirements.

## 2.7 Structural changes or variations

No arbitrary structural changes or variations to the Product are admitted. Any modifications must have the prior written authorization of TECNO-GAZ. In case of the Product having been tampered with, the warranty shall be invalidated and the manufacturer disclaims all liability for any injuries or damage caused to the OPERATOR, the RESPONSIBLE ORGANIZATION and the TECHNICAL SERVICE PERSONNEL.

## 3 Instructions on how to prepare the premises mechanically and electrically

### 3.1 Preparing the premises mechanically (ceiling and wall version)

The masonry works for preparing the ceiling to install the Product must be sturdy and safe and performed in a workmanlike manner according to applicable building regulations.

By way of example only, the professional persons charged with completing the masonry works are: Construction Engineer, Draughtsman, Building firm, duly registered in a professional register.

In case of wrong perforation of the Product supporting wall (e.g., the breakage of a reinforced-concrete ceiling iron) always inform the building manager as this could affect the stability of the building.

The ceiling must be able to withstand a weight of at least 300 kg/m<sup>2</sup> and have a thickness of at least 250 mm.

The Product installation premises must conform to local building standards.

After making sure the premises used for medical purposes are in conformity with the above requirements, proceed to mechanically anchor the ceiling plate, assessing the type of building and making all consequent adaptations.

The TECHNICAL SERVICE PERSONNEL have all technical, civil and legal responsibility relating to correctly and suitably performing Product anchoring and installation operations in a workmanlike manner.

## 3.2 Correctly wiring up the premises

The premises used for medical purposes must be safely wired up in a workmanlike manner by TECHNICAL SERVICE PERSONNEL to power the Product.

Before installing the Product, the TECHNICAL SERVICE PERSONNEL must make sure the following conditions exist:

- The wiring system of the environment (premises) in which installation is made must be in conformity with regulations for premises used for medical purposes and with applicable national laws and/or regulations.
- The electrical system must have a certificate of conformity issued by whosoever installed it.

The earth system must be certified as required by applicable regulations.

Position the thermal magnetic switch near the Product so that it can be switched off in case of need.

In case of mobile and wall versions do not position the device so it is hard to reach and remove the power plug in case of an emergency.


## 4 Product installation

**Before proceeding to install the Product, first of all check the presence of all the packaging and that this is in good condition and has not been damaged during transport.**

**Claims will only be taken into consideration if the seller or carrier has been immediately notified. All claims must be made in writing. Goods always travel under the responsibility and at the risk of the buyer.**

**Keep the original packaging in case the Product has to be re-dispatched.**

Personnel required:

Necessary protection equipment:

- Safety eyewear
- Gloves
- Accident-prevention footwear

Main switch


CAUTION

Ceiling version

Mobile version

Wall version

Fixing positions

Special equipment:

- Drill (ceiling and wall versions only)
- Set of hexagon spanners
- Screwdriver
- Circlip pliers
- Ladder (ceiling and wall versions only)
- Standard manual tools
- Set of drill bits (ceiling and wall versions only)

**After installation, the Product must be tested by Technical Service Personnel before being used.**

## 4.1 Parts included in the package

The Product is supplied complete with lamp head, sterilisable handpiece, swing arm, horizontal arm, bar, bar cover with relative safety ring, structure retention screws with glue, switchboard. TECNO-GAZ does not provide any kind of anchoring for fastening the plate to the ceiling. Such equipment must all be provided by the installer.

The Product is supplied complete with lamp head, sterilisable handpiece, swing arm, stems, wheeled base and base cover.


The Product is supplied complete with lamp head, sterilisable handpiece, swing arm, horizontal arm, wall plate with switchboard and HAM M6x50 wall plugs.

## 4.2 Ceiling and wall drilling instructions

For ceiling installation, the length of the bar varies according to the height of the premises in which the Product is installed.

The length of the bar is calculated to install the Product at a finished height off the floor of around 2000/2100 mm (as per drawing below), unless otherwise requested by the RESPONSIBLE ORGANIZATION.


Reinforced concrete  
Mechanical anchoring

By way of example only, below is a list of some types of walls:  
Proceed to fasten the ceiling plate using Hilti HSL-3-G M16/25 screw anchors (6 HAM M6x50 wall anchors) or other anchors with similar characteristics, carefully following the instructions provided by the anchor manufacturers and shown below for information purposes:

|    |                        |           |                |                        |
|----|------------------------|-----------|----------------|------------------------|
| do | Nominal diameter       | drill bit | Mt             | Closing bending moment |
| t  | Minimum drilling depth | Sw        | Wrench opening |                        |
| hs | Minimum depth          | insertion | x              | Fastening height       |
| l  | Anchor tie-rod length  |           |                |                        |

| Anchor tie-rod  | do<br>(mm) | t<br>(mm) | hs<br>(mm) | l<br>(mm) | Mt<br>(Nm) | SW<br>(mm) | x<br>(mm) |
|-----------------|------------|-----------|------------|-----------|------------|------------|-----------|
| HSL-3-G M 16/25 | 24         | 125       | 100        | 163       | 80         | 24         | 25        |
| HAM M6x50       | 12         | 65        | 40         | 50        | 10         | 10         | 20        |

1. Apply the paper template at the Product installation point and mark the fastening hole points with a pencil.
2. Make the holes in the ceiling in accordance with the anchoring tie-rod manufacturer's specifications.
3. Using a pump or a vacuum cleaner, remove the drilling residues and dust from the hole.
4. Fasten the Bar to the ceiling and, using a hammer, insert the anchoring tie-rod in the hole.  
**Attention!**  
**Check the fitting depth**
5. Using a torque wrench, tighten the anchorage to the tightening force indicated by the screw anchor manufacturer.  
The anchorage will immediately bear the weight.
6. Proceed in the same way for the remaining anchors.

## Chemical anchoring

## Hollow-core concrete

## Wall version

7. After one hour, again tighten the tie-rods to the prescribed tightening torque.

Drill the ceiling using the template provided. Insert the resin inside the 6 drilled holes and fill the hole, carefully following the manufacturer's instructions. TECNO-GAZ recommends using HILTI HIT-HY 270 resin or similar products.

Fit 6 threaded bars into the holes. TECNO-GAZ recommends M16 bars. Proceed to fasten the ceiling plate with nuts and locknuts for each tie-rod and tighten using the Allen key.

In this case the ceiling must be sandwich closed by means of the ceiling plate and counter-plate, being careful to include at least one rafter.

The plate and counter-plate must be fastened together using suitable M16 threaded steel bars, with ultimate tensile strength of at least 800 MPa, blocked at the top and bottom ends by relative washers, nuts and locknuts.

The Product is supplied complete with wall plate and HAM M6x50 wall plugs.


## 4.3 Instructions for ceiling version of Product

### 4.3.1 Installation of the ceiling plate, bar, power supply


VERSION WITH PRE-PREPARED TIE RODS OR CHEMICAL ANCHORAGE

If the tie rods are not prearranged in advance, place the template (drawing 12) (2) on the ceiling (1) and secure it with adhesive tape (3).

Drill the holes according to paragraph 4.2 and insert the 6 threaded bars M16 (4) into the ceiling.


**Make sure the product is stable.**

**Product falling hazard.**

On the threaded bars (or pre-prepared tie rods) insert the nuts (3), the washers (4) (secure them with adhesive tape on the tiges (bar) to prevent them from falling) and insert the tiges (2). Position the washers (5), nuts (6) and locknuts (7) from underneath. Using the nuts (6) and locknuts (7), position the tiges tube vertically, making sure it is correctly aligned using a spirit level (8). Once in the balanced position, tighten the upper nuts (3) so that everything is stable.


## VERSION WITH COUNTER-PLATE

In case of counter-plate (optional), drill the holes as shown in paragraph 4.2 and fix it to the ceiling.

To fasten the tiges to the counter-plate, follow the instructions given above.

Make sure the mains power cable (1) can reach the lamp power board without creating interferences with the Bar.

Fit the switchboard (2) on the Bar tube and tighten the two screws (3) and their toothed washers (4). Position the switchboard so the slot (5) of the retention bracket coincides with the hole M6 (6) on the fastening tube. Make sure the switchboard is secure by tightening the screw (7) and its toothed washer (8).

## 4.3.2 Installation of structure to bar

Align the pin of the horizontal arm (1) with the ceiling plate (2). Fit the connection cables (3) in the tube so that they come out from the side hole, in order to connect them to the power board. Insert the pin in the tube until the 3+3 holes at 120° of the pin coincide with the 3+3 holes at 120° of the tube.

- Insert all 6 screws (4).
- Strongly tighten ONLY two screws of a same side in vertical the one to the other.
- Now tighten the remaining screws.

This way the loosening over time will be avoided during continuous Product rotation.

**Product falling hazard.**

## 4.3.3 Installation of swinging arm



**INSERT THE PIVOT FROM THE BOTTOM  
TILL TO THE COMPLETE FIXING  
THROUGH THE 'QUICK COUPLING'**

### **WARNING, PERSONAL INJURY RISK**

**It's compulsory to follow this instruction.**

**The wrong insertion of locking levers in their seat of the pivot may cause the fall of swinging arm and cupola with high risk of personal injury.**


Align the swinging arm pin (1) and insert it into the horizontal arm (2) until the two locking levers (3) automatically engage and produce a "CLICK", and lock the arm in place. During insertion the assembly ring (4) will be automatically ejected. It is only designed to guide and facilitate assembly. It is recommended to keep it in case of future installation of the Product.

**If the assembly ring comes out before the pin is inserted, assembly is not possible. In the event of this coming out, it must be repositioned in its seat.**

**Before continuing with the assembly, check with swinging arm movements that the locking levers are in place in the pin slot.**

Join together connectors (1) and (2). In case of a standard lamp only one locking connector will be available. In case of a lamp equipped with CCTV, the supply will include power connectors, to be connected according to colours, and video signal connectors, to be connected according to letters. These connectors need to be screwed together. Put the wires into the horizontal arm slot (3). Place the plastic cap (4) on the upper part of the horizontal arm making sure that the 4 tips fit into their respective seats and close with the screws (5-6). Tighten the clutches (7) in order to make the arm position stable.

**Before proceeding, make sure the plastic cap is correctly positioned and in contact with the horizontal arm and that the screws are well tightened.**



**The correct closure of the plastic cap ensures that the locking levers are locked.**

**To make assembly easier, it is best to assemble the swinging arm first, and then the cupola.**

**Product falling hazard.**

## 4.3.4 Installation of cupola

Before positioning the cupola, as indicated in the drawing, position the yoke lock covering (1) and the yoke lock (2) on the swinging arm tube.

Insert the head of the Product with the yoke (3) on the swinging arm and bring the surfaces into contact. Now the head is able to maintain the position autonomously, without any support.

Pay attention to place the Product head and the arm in the same position, as indicated in the drawing, with the yoke to the left of the arm and the friction screw (4) turned downwards.


Push the yoke hub lock (2) onto the yoke hub (3) and rotate it in order to match the corresponding 6 holes.

Screw the 6 screws (5) to lock the hub and lock.


Then position the cover (1) on the lock (2) in order to cover the screws.

To ease the cables connection, first slightly extract the white connector coming from the swinging arm and then the one coming from the yoke.


Connect the white connectors together as shown in the images to the side


Then reposition the connectors inside the yoke, taking care not to crush cables.


In case of a lamp equipped with CCTV, in addition to the power connector, there will also be video signal connectors.

Connect the connectors in accordance with the letters.

Then also reposition these connectors inside the yoke, taking care not to crush cables.

Once the connection is complete, fasten the covering disc (1) in front of the yoke by screwing the three screws (2).

**Electric shock hazard.**

**For wiring connections in the ceiling version, use a cable suitable for at least 105°C and connect the ground lead to the terminal provided.**

**Strictly follow the wiring diagram for the correct connection.**

**Comply with the line and neutral connection indications. An inverted connection would lead to the lack of fuse protection.**

**Failure to connect the earth cables would prevent the safety of the Product.**

## 4.3.5 Electrical connection

Before completing the installation, it is necessary to connect the electrical system. To avoid any risk of electric shocks, the Product must only be connected to mains supplies with earth protection.

Before making the Product power connections, make sure the mains supply line has been interrupted.

TECNO-GAZ does not supply the mains supply cables.

Prime the fuses in the switchboard terminal box after the mechanical and electrical assembly of the Product. Priming the fuses too early could permanently damage the Product. If the Product is not used for long periods of time, remove the fuses.

According to the different cases, the electronic panel could be installed on the anchoring tube or in the false ceiling.

The line and neutral cables (L, N) from the power line have to be connected into the terminal 1 and 2 for single light head lamps and for double light heads to the terminals 1 and 2, 3 and 4.

Connect the wires (red and black) to the 3 and 4 terminals in case of single light head lamps and in case of double light head lamps to the terminals 5 and 6, 7 and 8. Follow always the colors and numbers on cables and terminals.

Always connect earthing cables (⏏) of lamp and net into the related terminals.


Out of every lamp there is always a communication cable (OW) connected to a forbox. In case of double light head lamp the cables (OW) of both lamps are connected to the same forbox.

This OW cable is needed for the communication between the lamps and to the optional wall control.

The cable (OW) of the wall control panel (if supplied) has to be connected to the forbox in order to allow the communication.

If the lamp isn't equipped with wall control panel, do not consider forbox connection.

### 4.3.6 Installation of ceiling cover

Once the electrical connections have been completed, the installation can be completed by positioning the ceiling cover. Depending on the type of ceiling (with false ceiling or not), a cover is provided split into two halves, which can be high or low. To install, position the two halves (1) in line with the ceiling anchoring tube. Close them by tightening the two screws of the ring (2) and the 4 screws of the cover (3).

Fasten the ring earth lead in the respective terminal.

Bring the complete cover (4) up against the ceiling / false ceiling and secure it in position by fully tightening the 4 screws (5).


## 4.4 Installation of Product in mobile version

### 4.4.1 Installation of lamp stem


Position the lower stem (2) in the base housing (1) and tighten it with the 4 screws (3).

Adequately tighten the 4 screws (3) to avoid any risk of instability and possible Product overturning.

Insert from the top of the stem (2) the stand cover (4), the closing ring (5) and the stem cover (6) in the indicated order.


**Instability and overturning hazard.**

Insert the cables inside the top stem (2) and position it vertically above the lower stem (1). Make the two extremities coincide using the guides (3).

Fasten the two stems by means of the screws (4).

## 4.4.2 Installation of swinging arm

Position the swinging arm (1) in front of the stem and in correspondence with the pivot (2).

Match the threaded hole of pivot (3) with the hole located on the hub (4).

Insert the swinging arm (1) into the pivot (2) and fastening it by tightening the screw (5).

Insert the plastic cover (6) from the bottom, widening it if required to make insertion easier. Fasten the cover by inserting the fasteners in the hub recesses.

Join the wiring connectors and fasten the upper cover (7) with the screw (8).

Position the cover (9) and secure it with the screw (10) in line with the threaded hole.

## 4.4.3 Installation of cupola


See the previous point 4.3.4.

## 4.4.4 Electrical connection

Lift the closing ring (1) and the stand cover (2) by 30-40 cm in order to access the power section. Join the connectors coming from the stem and switch. Return the cover and seal to original position and faster the cover (2) by means of the screws (3) to be fastened to the threaded bush (4). In the case of a battery lamp, also connect the battery faston that is disconnected.

After making the connection, engage the fuses.

The Product power connection is by means of a plug integrated in the supply cable supplied with the Product, placed on the box.


## 4.5 Installation of Product in wall version


### 4.5.1 Installation of plate with wall switchboard

Make sure the wall plate is fastened level, so the Product is in stable position.

Position the template sheet (2) (drawing 11) on the wall (1) and fasten it with adhesive tape (3), with the aid of a spirit level (4) to ensure levelling.

Make the holes as indicated in paragraph 3.1.


Fasten the plate (2) to the wall (1) with the aid of a spirit level (3).

**Product falling hazard.**

## 4.5.2 Installation of structure to plate

Insert in the wall box pin (1) a bronze washer (2), the Teflon washer (3) and the other bronze washer (2).

Afterwards, align the horizontal arm (4) with the box pin and insert it until the two locking levers (5) automatically engage and produce a "CLICK", and lock the arm in place.

During insertion the assembly ring (6) will be automatically ejected. It is only designed to guide and facilitate assembly.


It is recommended to keep it in case of future installation of the Product.

**If the assembly ring comes out before the pin is inserted, assembly is not possible. In the event of this coming out, it must be repositioned in its seat.**

**Before continuing with the assembly, check with swinging arm movements that the locking levers are in place in the pin slot.**

Join together connectors (1) and (2). In case of a standard lamp only one locking connector will be available.

In case of a lamp equipped with CCTV, the supply will include power connectors, to be connected according to colours, and video signal connectors, to be connected according to letters. These connectors need to be screwed together.

Put the wires into the horizontal arm slot (3).

Place the plastic cap (4) on the upper part of the horizontal arm making sure that the 4 tips fit into their respective seats and close with the screws (5-6).

Tighten the clutches (7) in order to make the arm position stable.

**Before proceeding, make sure the plastic cap is correctly positioned and in contact with the horizontal arm and that the screws are well tightened.**

**The correct closure of the plastic cap ensures that the locking levers are locked.**


**Product falling hazard.**

### 4.5.3 Installation of swinging arm

**INSERT THE PIVOT FROM THE BOTTOM  
TILL TO THE COMPLETE FIXING  
THROUGH THE 'QUICK COUPLING'**

**WARNING, PERSONAL INJURY RISK**  
It's compulsory to follow this instruction.  
The wrong inserction of locking levers in their seat of  
the pivot may cause the fall of swinging arm and  
cupola with high risk of personal injury.


Align the swinging arm pin (1) and insert it into the horizontal arm (2) until the two locking levers (3) automatically engage and produce a "CLICK", and lock the arm in place. During insertion the assembly ring (4) will be automatically ejected. It is only designed to guide and facilitate assembly. It is recommended to keep it in case of future installation of the Product.

**If the assembly ring comes out  
before the pin is inserted,  
assembly is not possible. In the  
event of this coming out, it must  
be repositioned in its seat.**

**Before continuing with the  
assembly, check with swinging  
arm movements that the locking  
levers are in place in the pin slot.**

**Before proceeding, make sure the plastic cap is correctly positioned and in contact with the horizontal arm and that the screws are well tightened.**

**The correct closure of the plastic cap ensures that the locking levers are locked.**

Join together connectors (1) and (2). In case of a standard lamp only one locking connector will be available.

In case of a lamp equipped with CCTV, the supply will include power connectors, to be connected according to colours, and video signal connectors, to be connected according to letters. These connectors need to be screwed together.

Put the wires into the horizontal arm slot (3).

Place the plastic cap (4) on the upper part of the horizontal arm making sure that the 4 tips fit into their respective seats and close with the screws (5-6).


Tighten the clutches (7) in order to make the arm position stable.

## 4.5.4 Installation of cupola

See the previous point 4.3.4.

## 4.5.5 Electrical connection

Remove the cover (1) of the wall box by loosening the screws (2). Engage the fuses and reposition the wall box cover (1) in place and tighten the screws (2).

## Fuses

**4.6 Protection fuses**

Product power protection is ensured by input fuses (L, N) and one output fuse (24V) of the TXAH 250V 5x20 type (where X is the fuse value). Only one input fuse (L) for ceiling versions.

FOR CEILING MODELS:

n°1 T2AH (L) and n°1 T10AH (+ 24VDC)

FOR MOBILE MODELS:

n°2 T2AH (L-N) and n°1 T10AH (+ 24VDC)

FOR MOBILE BATTERY MODEL:

n°2 T4AH (L-N) and n°1 T10AH (+ 24VDC)

FOR WALL MODEL:

n°1 T2AH (L-N) and n°1 T10AH (+ 24VDC)

**4.7 Handpiece fitting**

Insert the handpiece up tight on the support and turn it until the steel lever engages in its original place and rotation is blocked. Finally make sure the handpiece is well secured.

**4.8 Mechanical adjustments**

The Product is supplied correctly clutched and balanced. To make movement adjustment, refer to the setting instructions shown in the operation and maintenance manual.

**4.9 First switch-on**

To ensure the Product operates correctly, proceed as follows:

1. Make sure the power rating of the premises corresponds to that of the Product;
2. Fit the plug in the power socket of the premises – mobile and wall versions only;
3. Close the switch upstream of the system;
4. Move the Product switch located on the base cover for the mobile version and on the power box for the wall version respectively to position "I" (ON);
5. Press the  keyboard positioned on the lower part of the Product cupola.
6. Make sure all LEDs and functions are working properly.

At the time of commissioning, perform the electrical tests and prescriptions indicated in the IEC 62353 standard.

#### **4.10 Check the result of Product installation and testing before use**

Ticking the requirements listed below, if applicable to the Product version, is mandatory to ensure correct installation.

1. Make sure the ceiling is suitable for Product installation. ☐
2. Using a spirit level, make sure the bar is perpendicular with the ceiling or that the wall plate is horizontal with the wall. ☐
3. Make sure the switchboard is correctly fastened to the Bar by means of the threaded hole provided. ☐
4. Make sure the screws sustaining the horizontal arm are tight (*ceiling and wall versions*). ☐
5. Check that the locking levers are in place and the cap with the 4 tips is inserted correctly (*ceiling versions*). ☐
6. Make sure the stand is correctly fitted in the base (*mobile version*). ☐
7. Check the Product earth connection and make sure the earth terminals are well tightened. ☐
8. Check the correct rotation of the articulated joints and mechanical movements. ☐
9. Adjustment and rotation movements must be carefully clutched to ensure the Product is stable and maintains its position. ☐
10. Make sure the Product emits light. ☐

Stamp and signature of TECHNICAL SERVICE PERSONNEL:

-----

## 5 Troubleshooting

| N. | Problem                                               | Solution                                                                                                                                                                                    |
|----|-------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1  | The Product fails to remain in stable position        | Make sure the instructions in this manual, in the "Product installation" paragraph, have been correctly followed.<br>Refer to setting instructions in the operation and maintenance manual. |
| 2  | The Product fails to work                             | Make sure fuses have been fitted inside the terminal board. Make sure the electrical connectors are fitted. Check if there is voltage inside the Product.                                   |
| 3  | The fuse continues to burn out                        | Check the specifications of the fitted fuses.                                                                                                                                               |
| 4  | The light flickers and produces a stroboscopic effect | Contact the after sales service.                                                                                                                                                            |
| 5  | The Product does not switch on                        | Check the supply power voltage and check the fuses.<br>The electronics are faulty: contact the after sales service.                                                                         |

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